



V I S T E K

I N C O R P O R A T E D

Series 2500 Vibration Isolation Table - Size Selection Chart

Selecting the correct and complete Part Number is a three-step process.

First, determine the dimensions of the table you want and select the correct Part Number Prefix.

Part Number	Table Dimensions			*
	W	D	H	L, C, R
VIT2500-363030-*	36	30	30	C
VIT2500-363036-*	36	30	36	C
VIT2500-363630-*	36	36	30	C
VIT2500-363636-*	36	36	36	C
VIT2500-423030-*	42	30	30	L, C, R
VIT2500-423036-*	42	30	36	L, C, R
VIT2500-423630-*	42	36	30	L, C, R
VIT2500-423636-*	42	36	36	L, C, R
VIT2500-483030-*	48	30	30	L, C, R
VIT2500-483036-*	48	30	36	L, C, R
VIT2500-483630-*	48	36	30	L, C, R
VIT2500-483636-*	48	36	36	L, C, R
VIT2500-603030-*	60	30	30	L, C, R
VIT2500-603036-*	60	30	36	L, C, R
VIT2500-603630-*	60	36	30	L, C, R
VIT2500-603636-*	60	36	36	L, C, R
VIT2500-723030-*	72	30	30	L, C, R
VIT2500-723036-*	72	30	36	L, C, R
VIT2500-723630-*	72	36	30	L, C, R
VIT2500-723636-*	72	36	36	L, C, R

Second, determine whether the isolated area is to be in the center (C) of the table, on the left (L) hand side or right (R) hand side. The * in the Part Number designates where you designate your choice. The far right column above indicates which options are available for each table size.

Sample Part Numbers: VIT2500-363030-C
VIT2500-483630-L

Third, determine the weight of the instrument to be isolated, and select the corresponding part number below. If the weight of your instrument places you in two payload ranges, select the part number that places you towards the middle of a range.

Example: Your microscope weighs 45 Lbs. In this case, choose PN VIT2500-PL-02 and not VIT2500-PL-01.

This Part Number does not represent a product. It simply tells us what bearings to use on the table selected with the Part Number from the chart in the left hand column of this page.

Payload Part Number	Instrument Weight Range	
	Min. Lbs (Kg)	Max. Lbs (Kg)
VIT2500-PL-01	5 (2)	60 (27)
VIT2500-PL-02	30 (14)	70 (32)
VIT2500-PL-03	50 (23)	105 (48)
VIT2500-PL-04	65 (29)	125 (57)
VIT2500-PL-05	95 (43)	150 (68)
VIT2500-PL-06	120 (54)	180 (82)
VIT2500-PL-07	145 (66)	275 (125)

Vistek Inc.
1741 W, University Drive, Ste 146.
Tempe, AZ 85281
P.602.426.1211/F.602.426.1611
www.vistekinc.com
